

Package ‘quicR’

February 26, 2025

Title RT-QuIC Data Formatting and Analysis

Version 2.1.0

Description Designed for the curation and analysis of data generated from real-time quaking-induced conversion (RT-QuIC) assays first described by Atarashi et al. (2011) <[doi:10.1038/nm.2294](https://doi.org/10.1038/nm.2294)>. ‘quicR’ calculates useful metrics such as maxpoint ratio: Rowden et al. (2023) <[doi:10.1099/vir.0.069906-0](https://doi.org/10.1099/vir.0.069906-0)>; time-to-threshold: Shi et al. (2013) <[doi:10.1186/2051-5960-1-44](https://doi.org/10.1186/2051-5960-1-44)>; and maximum slope. Integration with the output from plate readers allows for seamless input of raw data into the R environment.

Imports dplyr, ggplot2, janitor, openxlsx, purrr, readxl, reshape2, slider, stats, stringr, tidyverse

Depends R (>= 4.1.0)

License GPL-3

Encoding UTF-8

RoxygenNote 7.3.2

Suggests testthat (>= 3.1.0)

Config/testthat/edition 3

NeedsCompilation no

Author Gage Rowden [aut, cre] (<<https://orcid.org/0000-0002-7517-0480>>)

Maintainer Gage Rowden <gage.rowden1145@gmail.com>

Repository CRAN

Date/Publication 2025-02-26 19:30:02 UTC

Contents

| | |
|-------------------|---|
| quicR-package | 2 |
| add_reps | 2 |
| BMG_format | 3 |
| calculate_metrics | 4 |
| calculate_MPR | 5 |
| calculate_MS | 6 |

| | |
|--------------------------------|----|
| calculate_threshold | 6 |
| calculate_TtT | 7 |
| convert_tables | 8 |
| get_meta | 9 |
| get_real | 9 |
| get_sample_locations | 10 |
| get_wells | 11 |
| normalize_RFU | 12 |
| organize_tables | 13 |
| plate_view | 13 |
| plot_metrics | 14 |
| separate_raw | 15 |
| transpose_real | 16 |

| | |
|---------------|---|
| quicR-package | <i>@description Data handling for real-time quaking induced conversion assays</i> |
|---------------|---|

Description

See the README on [Github](#)

Details

@keywords internal

Author(s)

Maintainer: Gage Rowden <gage.rowden1145@gmail.com> ([ORCID](#))

| | |
|----------|-----------------------|
| add_reps | <i>Add replicates</i> |
|----------|-----------------------|

Description

Adds replicate information to the sample IDs. Well IDs should be formatted like so: A4, B9, H11, J24

Usage

```
add_reps(df, sep = "_")
```

Arguments

- df A dataframe containing two columns for well IDs and Sample IDs
 sep a character string to separate the terms.

Value

A dataframe with replicate numbers pasted to the Sample IDs

BMG_format

*Format Table for BMG Sample ID Import***Description**

BMG_format accepts a plate layout .CSV file and formats the Sample IDs into a format which can be easily imported into the BMG control software.

Usage

```
BMG_format(  
  file,  
  save_path = "",  
  save_name = "formatted.txt",  
  write_file = FALSE  
)
```

Arguments

- file A .CSV file containing the plate layout of Sample IDs.
 save_path The path to the directory that you want the file saved.
 save_name The name of the output file. Should have the ".txt" extension.
 write_file Logical. If true, function will write a .txt file; otherwise it will return a character vector.

Value

A text file containing information for import into the BMG control software.

Examples

```
layout_file <- system.file(  
  "extdata/BMG_formatting",  
  file = "plate_layout.csv",  
  package = "quicR"  
)  
BMG_format(layout_file)
```

| | |
|--------------------------------|--|
| <code>calculate_metrics</code> | <i>Generate a dataframe with calculated metrics.</i> |
|--------------------------------|--|

Description

Uses functions from the "calculate" family of quicR functions to generate an analyzed dataframe.

Usage

```
calculate_metrics(
  data,
  meta,
  metrics = c("MPR", "MS", "TtT", "RAF"),
  transpose = FALSE,
  normalize = FALSE,
  start_col = 3L,
  MS_window = 3L,
  threshold = 2
)
```

Arguments

| | |
|------------------------|--|
| <code>data</code> | A dataframe containing the raw RT-QuIC data. |
| <code>meta</code> | A dataframe containing sample metadata. Should include at least the "Sample IDs" column. |
| <code>metrics</code> | An array containing the metrics which should be calculated. |
| <code>transpose</code> | Logical; should the raw data be transposed before performing the calculations? |
| <code>normalize</code> | Logical; should the raw data be normalized before performing the calculations? |
| <code>start_col</code> | Integer; column number denoting where the numeric data begins. |
| <code>MS_window</code> | Integer; width of the window applied in the calculation of max slope. |
| <code>threshold</code> | Float; the threshold applied to the calculation of time-to-threshold. |

Value

A dataframe of calculated metrics.

Examples

```
file <- system.file(
  "extdata/input_files",
  file = "test4.xlsx",
  package = "quicR"
)

data <- quicR::get_real(file)[[1]] |>
  quicR::normalize_RFU()
```

```
meta <- quicR::organize_tables(file) |>  
quicR::convert_tables()  
  
calculate_metrics(data, meta)
```

| | |
|---------------|-------------------------------------|
| calculate_MPR | <i>Calculate the Maxpoint Ratio</i> |
|---------------|-------------------------------------|

Description

Maxpoint ratio is defined as the maximum relative fluorescence divided by the background fluorescence.

Usage

```
calculate_MPR(data, start_col = 3, data_is_norm = TRUE)
```

Arguments

| | |
|--------------|--|
| data | A dataframe containing the real-time fluorescence data. |
| start_col | Integer, the column at which the background fluorescence should be read. |
| data_is_norm | Logical, if the data has not been normalized, will make a call to normalize_RFU. |

Value

A vector containing MPR values.

Examples

```
# This test takes >5 sec  
  
file <- system.file(  
  "extdata/input_files",  
  file = "test.xlsx",  
  package = "quicR"  
)  
df_ <- quicR::get_real(file)[[1]]  
print(calculate_MPR(df_))
```

| | |
|---------------------------|--------------------------------|
| <code>calculate_MS</code> | <i>Calculate Maximum Slope</i> |
|---------------------------|--------------------------------|

Description

Uses a sliding window to calculate the slope of real-time reads.

Usage

```
calculate_MS(data, window = 3, data_is_norm = TRUE)
```

Arguments

| | |
|---------------------------|---|
| <code>data</code> | A dataframe containing real-time reads. It is recommended to use a dataframe made from <code>normalize_RFU</code> . |
| <code>window</code> | Integer designating how wide you want the sliding window to be for calculating the moving average slope. |
| <code>data_is_norm</code> | Logical; if FALSE, will make a call to <code>normalize_RFU</code> . |

Value

A dataframe containing the real-time slope values as change in RFU/sec.

Examples

```
# This test takes >5 sec

file <- system.file(
  "extdata/input_files",
  file = "rt_data.csv",
  package = "quicR"
)
df_ <- read.csv(file, check.names = FALSE)
calculate_MS(df_)
```

| | |
|----------------------------------|---|
| <code>calculate_threshold</code> | <i>Calculate a Threshold for Rate Determination</i> |
|----------------------------------|---|

Description

Calculates a threshold for determining time-to-threshold and rate of amyloid formation.

Usage

```
calculate_threshold(
  data,
  background_cycle = 2,
  method = list("stdev", "none"),
  multiplier = 1
)
```

Arguments

| | |
|------------------|---|
| data | A dataframe output from get_real. |
| background_cycle | Integer; the cycle used for background fluorescence. |
| method | Method for determining threshold; default is "stdev". |
| multiplier | For some methods, will add a multiplier for more conservative thresholds. |

Value

A float value.

Examples

```
file <- system.file(
  "extdata/input_files",
  file = "test2.xlsx",
  package = "quicR"
)
threshold <- get_real(file)[[1]] |>
  calculate_threshold(multiplier = 10)
```

calculate_TtT

Calculate Time to Threshold

Description

Calculates the time required to reach a defined threshold.

Usage

```
calculate_TtT(data, threshold, start_col = 3)
```

Arguments

| | |
|-----------|--|
| data | A dataframe containing real-time RT-QuIC data. |
| threshold | A numeric value defining the threshold. |
| start_col | The column containing the starting position of the real-time data. |

Value

A vector containing the times to threshold

Examples

```
# This test takes >5 sec

file <- system.file(
  "extdata/input_files",
  file = "test2.xlsx",
  package = "quicR"
)
df_ <- get_real(file)[[1]] |>
  quicR::transpose_real() |>
  quicR::normalize_RFU(transposed = TRUE)
calculate_TtT(df_, threshold = 2)
```

convert_tables

Convert tables into a single column in a dataframe.

Description

Accepts a table or matrix or a list of tables and matrices and converts them into dataframe columns.

Usage

```
convert_tables(tab, na.omit = TRUE)
```

Arguments

| | |
|---------|--|
| tab | A table/matrix or a list of tables/matrices. |
| na.omit | Logical; if true, will remove rows with NA. |

Value

A dataframe column.

Examples

```
file <- system.file(
  "extdata/input_files",
  file = "test.xlsx",
  package = "quicR"
)
tabs <- organize_tables(file)
convert_tables(tabs)
```

| | |
|----------|----------------------------------|
| get_meta | <i>Retrieve the BMG metadata</i> |
|----------|----------------------------------|

Description

Takes the Excel file exported from MARS and compiles the metadata in the header.

Usage

```
get_meta(file)
```

Arguments

file The Excel file exported from MARS.

Value

A dataframe containing the Meta_ID and Meta_info

Examples

```
file <- system.file(  
  "extdata/input_files",  
  file = "test.xlsx",  
  package = "quicR"  
)  
get_meta(file)
```

| | |
|----------|--|
| get_real | <i>Get Real-Time RT-QuIC Fluorescence Data</i> |
|----------|--|

Description

Accepts an Excel file or a dataframe of real-time RT-QuIC data.

Usage

```
get_real(data, ordered = FALSE)
```

Arguments

data Either an Excel file or a dataframe.

ordered Logical, if true, will organize the columns by sample ID rather than by well.

Value

A list of dataframes containing the formatted real-time data.

Examples

```
file <- system.file(
  "extdata/input_files",
  file = "test.xlsx",
  package = "quicR"
)
get_real(file)
```

`get_sample_locations` *Get the well locations of the samples used in the RT-QuIC run.*

Description

Returns a dataframe with the sample IDs and well IDs used in the plate.

Usage

```
get_sample_locations(
  file,
  tab_name = "Sample IDs",
  dilution_bool = FALSE,
  dilution_fun = function(x) 1 * x,
  sep = "\n",
  plate = 96
)
```

Arguments

| | |
|----------------------------|---|
| <code>file</code> | Excel file exported from MARS |
| <code>tab_name</code> | Table name containing the sample IDs. |
| <code>dilution_bool</code> | Logical; is there a table containing dilution factors? If so, will add a newline and the log of the dilution factor to the ID column. |
| <code>dilution_fun</code> | A function for transforming the dilution factor. |
| <code>sep</code> | A string used to separate the sample ID and dilution factor. |
| <code>plate</code> | Integer; either 96 or 384 to denote microplate type. |

Value

A vector containing well IDs.

Examples

```
file <- system.file(  
  "extdata/input_files",  
  file = "test.xlsx",  
  package = "quicR"  
)  
get_sample_locations(file)
```

get_wells*Get the Wells Used in the RT-QuIC Run.*

Description

Returns the well IDs used in the plate.

Usage

```
get_wells(file)
```

Arguments

file Excel file exported from MARS

Value

A vector containing well IDs.

Examples

```
file <- system.file(  
  "extdata/input_files",  
  file = "test.xlsx",  
  package = "quicR"  
)  
get_wells(file)
```

| | |
|---------------|-------------------------------|
| normalize_RFU | <i>Normalize Fluorescence</i> |
|---------------|-------------------------------|

Description

Normalizes the real-time RT-QuIC data against the background fluorescence of a defined cycle. All cycles are divided by the fluorescent value of the defined cycle.

Usage

```
normalize_RFU(data, bg_cycle = 4, transposed = FALSE)
```

Arguments

| | |
|------------|--|
| data | A dataframe generated from get_real. |
| bg_cycle | The cycle used for background fluorescence |
| transposed | Logical, TRUE if cycle values are shown as column names. |

Value

A dataframe containing real-time normalized fluorescence values.

Examples

```
# This test takes >5 sec

file <- system.file(
  "extdata/input_files",
  file = "test2.xlsx",
  package = "quicR"
)
df_ <- get_real(file)[[1]]

# Export the tables in the first sheet of the file.
dic <- quicR::organize_tables(file)

# Normalize the raw data against the background reading.
normalize_RFU(df_)
```

| | |
|-----------------|-----------------------------|
| organize_tables | <i>Organize MARS Tables</i> |
|-----------------|-----------------------------|

Description

Extracts the tables from the microplate view sheet in the MARS Excel file and adds each table to a list.

Usage

```
organize_tables(file, plate = 96)
```

Arguments

| | |
|-------|---|
| file | An Excel file exported from MARS. |
| plate | Integer either 96 or 384 to denote microplate type. |

Value

A list containing tibbles.

Examples

```
file <- system.file(  
  "extdata/input_files",  
  file = "test.xlsx",  
  package = "quicR"  
)  
organize_tables(file)
```

| | |
|------------|-----------------------------|
| plate_view | <i>Real-Time Plate View</i> |
|------------|-----------------------------|

Description

Converts the real-time data into a ggplot figure. The layout is either 8x12 or 16x24 for 96- and 384-well plates, respectively.

Usage

```
plate_view(df, meta, plate = 96)
```

Arguments

| | |
|--------------------|---|
| <code>df</code> | Real-time dataframe |
| <code>meta</code> | Dataframe containing well IDs and Sample IDs to title each facet. |
| <code>plate</code> | Integer either 96 or 384 to denote microplate type. |

Value

A ggplot object

Examples

```
# This test takes >5 sec

file <- system.file(
  "extdata/input_files",
  file = "test2.xlsx",
  package = "quicR"
)

# Get the real-time data.
df_ <- get_real(file, ordered = FALSE)[[1]] |>
  as.data.frame()

sample_locations <- get_sample_locations(
  file,
  dilution_bool = TRUE,
  dilution_fun = function(x) -log10(x)
)

plate_view(df_, sample_locations)
```

plot_metrics

Plot metrics generated from the "calculate" family of quicR functions.

Description

Generates a faceted figure of boxplots.

Usage

```
plot_metrics(
  data,
  sample_col = "Sample IDs",
  fill = "Dilutions",
  dilution_bool = TRUE,
  nrow = 2,
  ncol = 2
)
```

Arguments

| | |
|----------------------------|---|
| <code>data</code> | A dataframe containing the calculated metrics from the "calculate" family of quicR functions. |
| <code>sample_col</code> | The name of the column containing the sample IDs. |
| <code>fill</code> | The column containing the fill aesthetic. Usually the dilutions column. |
| <code>dilution_bool</code> | Logical; should dilution factors be included in the plot? |
| <code>nrow</code> | Integer; number of rows to output in the plot. |
| <code>ncol</code> | Integer; number of columns to output in the plot. |

Value

A ggplot object

Examples

```
# This test takes >5 sec

file <- system.file(
  "extdata/input_files",
  file = "test4.xlsx",
  package = "quicR"
)

data <- quicR::get_real(file)[[1]] |>
  quicR::normalize_RFU()

meta <- quicR::organize_tables(file) |>
  quicR::convert_tables()

calculate_metrics(data, meta) |>
  plot_metrics()
```

`separate_raw`

Separate Real-Time Data into separate dataframes.

Description

If multiple real-time reads were exported from MARS, `separate_raw` will parse them out and separate them. It will also export to an Excel file with each real-time data having its own sheet.

Usage

```
separate_raw(file, num_rows, export_name)
```

Arguments

- `file` An Excel file exported from MARS.
`num_rows` Number of rows in the header to ignore.
`export_name` The name of the original file or an original name.

Value

An Excel file with separated raw real-time data.

`transpose_real` *Transpose Real-Time Data*

Description

Transposes the real-time data table exported by the BMG software. Accepts output from the function, "get_real".

Usage

```
transpose_real(data)
```

Arguments

- `data` A dataframe generated from get_real.

Value

A transposed dataframe containing real-time normalized fluorescence values.

Index

add_reps, 2
BMG_format, 3
calculate_metrics, 4
calculate_MPR, 5
calculate_MS, 6
calculate_threshold, 6
calculate_TtT, 7
convert_tables, 8
get_meta, 9
get_real, 9
get_sample_locations, 10
get_wells, 11
normalize_RFU, 12
organize_tables, 13
plate_view, 13
plot_metrics, 14
quicR(quicR-package), 2
quicR-package, 2
separate_raw, 15
transpose_real, 16